

Product Datasheet

Universal SYBR Green qPCR Master Mix (orb2946035)

Catalog Number	orb2946035
Category	Tools
Description	This 2X universal SYBR Green qPCR Master Mix is a ready-to-use solution containing a hot-start DNA polymerase, SYBR Green I dye, dNTPs, Mg ²⁺ , and ROX reference dye. It enables sensitive and reproducible gene expression quantification and genotyping in diverse in vitro research applications, including molecular biology and diagnostics, by simplifying reaction setup and minimizing contamination risks.
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)