

Product Datasheet

ALG.APV-527 Antibody (orb2949950)

Catalog Number	orb2949950
Category	Antibodies
Description	ALG.APV-527 (APV-527) is a bispecific antibody targeting 4-1BB (CD137) and 5T4. It is designed to activate T cells and natural killer (NK) cells through 4-1BB only when it concurrently binds to tumor cells expressing 5T4. ALG.APV-527 demonstrates potent antitumor activity.
Target	Others
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)