

Product Datasheet

HX009 Antibody (orb2949973)

Catalog Number	orb2949973
Category	Antibodies
Description	HX009 is a bispecific antibody with attenuated CD47 binding that blocks both PD1 and CD47 pathways. It is a research tool for studying non-Hodgkin lymphoma (NHL), applicable in both in vitro cellular assays and in vivo preclinical models to investigate dual-targeting immunotherapy.
Target	Immunology/Inflammation related,PD-1/PD-L1
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	34.25 mg/mL
Purity	95.00% (May vary between batches)
MW	173.16 kDa
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)