

Product Datasheet

Rilvegostomig Antibody (orb2949979)

Catalog Number	orb2949979
Category	Antibodies
Description	Rilvegostomig is a humanized bispecific IgG1 antibody that simultaneously targets the PD-1 and TIGIT immune checkpoint receptors. It is designed to enhance T-cell activation and is investigated in preclinical and clinical research, particularly for applications in oncology such as non-small cell lung cancer (NSCLC).
Target	PD-1/PD-L1
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2640305-01-5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)