

Product Datasheet

RO7122290 Antibody (orb2949981)

Catalog Number	orb2949981
Category	Antibodies
Description	RO7122290 is a bispecific FAP-4-1BB-targeting fusion protein with immunomodulatory and antitumor potential for advanced solid tumors. Its application in preclinical research includes investigating T-cell activation in vitro and antitumor efficacy in murine models, supporting its development as an immunotherapeutic agent.
Target	Others
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)