

Product Datasheet

Forimtamig Antibody (orb2950013)

Catalog Number	orb2950013
Category	Antibodies
Description	Forimtamig (RG-6324) is a bispecific T-cell engager with a 2+1 format, targeting GPRC5D and CD3. Its Fc domain contains P329G LALA mutations to minimize effector function. This construct is utilized in preclinical research for investigating multiple myeloma, including in vitro cytotoxicity assays and in vivo xenograft models.
Target	Integrin
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)