

Product Datasheet

PIR Rabbit Polyclonal Antibody (orb2951427)

Catalog Number	orb2951427
Category	Antibodies
Description	PIR Rabbit Polyclonal Antibody is an unconjugated polyclonal antibody. It is reactive with Human samples and is suitable for use in ELISA, IHC and WB applications. The antibody was generated using an E. coli - derived recombinant Human PIR (Met1-Asn290) as an immunogen. It is supplied as a liquid and purified by antigen affinity column.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Immunogen	E. coli - derived recombinant Human PIR (Met1-Asn290).
UniProt ID	O00625
Tested applications	ELISA, IHC, WB
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Western blot analysis was performed using anti-PIR polyclonal antibody at 1 µg/mL on various samples. Lane 1: HepG2 cell lysate, Lane 2: K562 cell lysate, Lane 3: NIH3T3 cell lysate, Lane 4: A431 cell lysate, Lane 5: A375 cell lysate, Lane 6: MCF-7 cell lysate, Lane 7: Mouse liver lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)