

## Product Datasheet

### UXT Rabbit Polyclonal Antibody (orb2953038)

|                             |  |
|-----------------------------|--|
| <b>Catalog Number</b>       | orb2953038   |
| <b>Category</b>             | Antibodies   |
| <b>Description</b>          | Anti-UXT Polyclonal Antibody is a rabbit polyclonal antibody detecting ART-27 in ELISA, IHC, WB. Suitable for Human, Mouse, and Rat. Highlights Affinity Purified: Minimal background and high purity for reliable results. Multi-Application: Validated across multiple applications. Multi-Species: Cross-reactive for translational research. |
| <b>Clonality</b>            | Polyclonal   |
| <b>Species/Host</b>         | Rabbit   |
| <b>Isotype</b>              | IgG  |
| <b>Conjugation</b>          | Unconjugated   |
| <b>Reactivity</b>           | Human, Mouse, Rat  |
| <b>Form/Appearance</b>      | Liquid   |
| <b>Concentration</b>        | Lot dependent  |
| <b>Buffer/Preservatives</b> | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.  |
| <b>Purification</b>         | Purified by antigen affinity column.   |
| <b>Immunogen</b>            | E. coli - derived recombinant Human UXT (Glu15-His157).  |
| <b>UniProt ID</b>           | <b>Q9UBK9</b>  |
| <b>Tested applications</b>  | ELISA, IHC, WB   |
| <b>Antibody Type</b>        | Primary Antibody   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |   |
|------------------------|---|
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |



Western blot analysis was performed using anti-UXT polyclonal antibody at 1  $\mu$ g/mL on recombinant UXT (Catalog No. orb2967196).

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558