

Product Datasheet

DLAT Rabbit Polyclonal Antibody (orb2953132)

Catalog Number	orb2953132
Category	Antibodies
Description	Anti-DLAT Polyclonal Antibody is a rabbit polyclonal antibody detecting DLTA in ELISA, IHC, WB. Suitable for Human, Mouse, and Rat. Highlights Affinity Purified: Minimal background and high purity for reliable results. Multi-Application: Validated across multiple applications. Multi-Species: Cross-reactive for translational research.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Immunogen	E. coli - derived recombinant Human DLAT (Gln374-Trp631).
UniProt ID	P10515
Tested applications	ELISA, IHC, WB
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Western blot analysis was performed using anti-DLAT polyclonal antibody at 1 µg/mL on various samples. Lane 1: Hela cell lysate, Lane 2: MCF-7 cell lysate, Lane 3: Lncap cell lysate, Lane 4: K562 cell lysate, Lane 5: HCT116 cell lysate, Lane 6: NIH3T3 cell lysate, Lane 7: Mouse brain lysate, Lane 8: Mouse liver lysate, Lane 9: Mouse kidney lysate, Lane 10: Rat brain lysate, Lane 11: Rat liver lysate, Lane 12: Rat kidney lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)