

Product Datasheet

SIGIRR Rabbit Polyclonal Antibody (orb2954798)

| | |
|-----------------------------|---|
| Catalog Number | orb2954798 |
| Category | Antibodies |
| Description | Anti-SIGIRR Polyclonal Antibody is a rabbit polyclonal antibody detecting TIR8 in ELISA, IHC, WB. Suitable for Human. Highlights Affinity Purified: Minimal background and high purity for reliable results. Multi-Application: Validated across multiple applications. |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Concentration | Lot dependent |
| Buffer/Preservatives | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300. |
| Purification | Purified by antigen affinity column. |
| Immunogen | E. coli - derived recombinant Human SIGIRR (Met1-His118). |
| UniProt ID | Q6IA17 |
| Tested applications | ELISA, IHC, WB |
| Antibody Type | Primary Antibody |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Expiration Date

12 months from date of receipt.



Western blot analysis was performed using anti-SIGIRR polyclonal antibody at 1 µg/mL on recombinant SIGIRR (Catalog No. orb2967724).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558