

Product Datasheet

AIF1/IBA1 Rabbit Polyclonal Antibody (orb2955075)

| | |
|-----------------------------|--|
| Catalog Number | orb2955075 |
| Category | Antibodies |
| Description | Anti-AIF1/IBA1 Polyclonal Antibody is a rabbit polyclonal antibody detecting Ionized calcium-binding adapter molecule 1 in ELISA, IHC, WB, IF. Suitable for Mouse, Bovine, Human, Macaca mulatta, and Rat. Highlights Affinity Purified: Minimal background and high purity for reliable results. Multi-Application: Validated across multiple applications. Multi-Species: Cross-reactive for translational research. |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | Unconjugated |
| Reactivity | Bovine, Human, Monkey, Mouse, Porcine, Rat |
| Form/Appearance | Liquid |
| Concentration | Lot dependent |
| Buffer/Preservatives | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300. |
| Purification | Purified by antigen affinity column. |
| Immunogen | E. coli - derived recombinant Human AIF1/IBA1 (Met1-Pro147). |
| UniProt ID | P55008 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|----------------------------|---|
| Tested applications | ELISA, IF, IHC, WB |
| Antibody Type | Primary Antibody |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis was performed using anti-AIF1 polyclonal antibody at 1 µg/mL on various samples. Lane 1: THP-1 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)