

Product Datasheet

CD31/PECAM1 Monoclonal Mouse Monoclonal Antibody (orb2956219)

Catalog Number	orb2956219
Category	Antibodies
Description	Anti-Human CD31/PECAM1 Monoclonal Antibody (B3C380) is a mouse monoclonal antibody detecting PECAM-1 in ELISA, IHC, WB. Suitable for Human. Highlights Recombinant: Superior batch-to-batch reproducibility with defined sequences. Multi-Application: Validated across multiple applications.
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	Liquid in PBS containing 50% Glycerol, 0.5% rAlbumin and 0.02% Sodium azide.
Purification	Protein A/G purified from cell culture supernatant.
Immunogen	Mammalian cells - derived recombinant Human CD31 / PECAM1 (Met1-Lys601).
UniProt ID	P16284
Tested applications	ELISA, IHC, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Clone Number	B3C380
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



SDS-PAGE for Anti-Human CD31/PECAM1 Monoclonal Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)