

Product Datasheet

CD284/TLR4 Monoclonal Mouse Monoclonal Antibody (orb2956314)

Catalog Number	orb2956314
Category	Antibodies
Description	Anti-Human CD284/TLR4 Monoclonal Antibody (B3C279) is a mouse monoclonal antibody detecting Toll-like receptor 4 in ELISA, IHC, WB. Suitable for Human. Highlights Recombinant: Superior batch-to-batch reproducibility with defined sequences. Multi-Application: Validated across multiple applications.
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG3
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Immunogen	Mammalian cells - derived recombinant Human CD284 / TLR4 (Met1-Lys631).
Tested applications	ELISA, IHC, WB
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	B3C279
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



SDS-PAGE for Anti-Human CD284/TLR4 Monoclonal Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)