

## Product Datasheet

### Ara h 2 Human IgE Antibody (38B7) (orb2956705)

|                             |                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2956705                                                                                                                              |
| <b>Category</b>             | Antibodies                                                                                                                              |
| <b>Description</b>          | Ara h 2 Human IgE Antibody (38B7)                                                                                                       |
| <b>Clonality</b>            | Monoclonal                                                                                                                              |
| <b>Species/Host</b>         | Human                                                                                                                                   |
| <b>Isotype</b>              | IgE                                                                                                                                     |
| <b>Conjugation</b>          | Unconjugated                                                                                                                            |
| <b>Reactivity</b>           | Plant                                                                                                                                   |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                  |
| <b>Concentration</b>        | 1 mg/ml                                                                                                                                 |
| <b>Buffer/Preservatives</b> | 0.01M PBS, pH 7.4.                                                                                                                      |
| <b>Purity</b>               | >95% as determined by SDS-PAGE.                                                                                                         |
| <b>Tested applications</b>  | ELISA                                                                                                                                   |
| <b>Antibody Type</b>        | Recombinant Antibody                                                                                                                    |
| <b>Endotoxins</b>           | Please contact with the lab for this information.                                                                                       |
| <b>Storage</b>              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                 | For research use only                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)