

Product Datasheet

HCMV/HHV5 gB Neutralizing Human Monoclonal Antibody (orb2956879)

Catalog Number	orb2956879
Category	Antibodies
Description	HCMV/HHV5 gB Neutralizing Human Monoclonal Antibody is an unconjugated monoclonal antibody. It shows reactivity with Virus samples. It is suitable for use in ELISA, NeA applications. The antibody is supplied in liquid form. It is Protein A/G purified from cell culture supernatant. Purity is >95% as determined by SDS-PAGE.
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG1, kappa
Conjugation	Unconjugated
Reactivity	Virus
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
UniProt ID	P13201
Tested applications	ELISA, NeA

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Antibody Type	Recombinant Antibody
Clone Number	BCD5X47
Endotoxins	Please contact with the lab for this information.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



SDS-PAGE for InVivoMAb Anti-HCMV/HHV5 gB Neutralizing Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558