

## Product Datasheet

### PD-1/PD-L1 Mouse Monoclonal Antibody (orb2956958)

|                      |                                                                                                                                                   |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number       | orb2956958                                                                                                                                        |
| Category             | Antibodies                                                                                                                                        |
| Description          | InVivoMAb Anti-Mouse PD-1 & PD-L1 Bispecific Antibody (BKX2462) is a mouse monoclonal antibody detecting Bispecific in FuncS. Suitable for Mouse. |
| Clonality            | Recombinant                                                                                                                                       |
| Species/Host         | Mouse                                                                                                                                             |
| Isotype              | IgG                                                                                                                                               |
| Conjugation          | Unconjugated                                                                                                                                      |
| Reactivity           | Mouse                                                                                                                                             |
| Form/Appearance      | Liquid                                                                                                                                            |
| Concentration        | Lot dependent                                                                                                                                     |
| Buffer/Preservatives | 0.01M PBS, pH 7.4.                                                                                                                                |
| Purity               | >95% as determined by SDS-PAGE.                                                                                                                   |
| Purification         | Protein A/G purified from cell culture supernatant.                                                                                               |
| UniProt ID           | <b>Q02242, Q9EP73</b>                                                                                                                             |
| Tested applications  | Functional Studies                                                                                                                                |
| Antibody Type        | Recombinant Antibody                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom


Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |                                                                                                                                         |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clone Number</b>    | BKX2462                                                                                                                                 |
| <b>Endotoxins</b>      | Please contact with the lab for this information.                                                                                       |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only                                                                                                                   |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                         |



SDS-PAGE for InVivoMAb Anti-Mouse PD-1 & PD-L1 Bispecific Antibody.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)