

Product Datasheet

CD79b/CD3 Monoclonal Antibody (orb2957212)

Catalog Number	orb2957212
Category	Antibodies
Description	InVivoMAb Anti-Human CD79b & CD3 Bispecific Antibody (BDU223) is a human monoclonal antibody detecting B-cell-specific glycoprotein B29 in FuncS. Suitable for Human.
Clonality	Recombinant
Species/Host	Human
Isotype	IgG like Duobody
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
UniProt ID	P40259, P07766
Tested applications	Functional Studies
Antibody Type	Recombinant Antibody
Clone Number	BDU223

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Endotoxins

Please contact with the lab for this information.

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only

Expiration Date

12 months from date of receipt.

**Anti-Human CD79b & CD3 Bispecific Antibody
binds with CD79b**

Detects recombinant CD79b (Catalog No: orb2968587) in indirect ELISA.




SDS-PAGE for InVivoMAb Anti-Human CD79b & CD3 Bispecific Antibody.

Biorbyt Ltd.7 Signet Court, Swann Road
Cambridge
CB5 8LA
United KingdomEmail: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United StatesEmail: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Western blot analysis was performed using anti-CD3E monoclonal antibody at 1 $\mu\text{g}/\text{mL}$ on various samples. Lane 1: recombinant human CD3E.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)