

Product Datasheet

GAD65/67 Mouse Monoclonal Antibody (orb2958060)

| | |
|-----------------------------|---|
| Catalog Number | orb2958060 |
| Category | Antibodies |
| Description | Anti-GAD65/67 Antibody (BN349X4) is a mouse monoclonal antibody detecting Glutamate decarboxylase 2 in IHC, WB. |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Isotype | IgG2a, kappa |
| Conjugation | Unconjugated |
| Form/Appearance | Liquid |
| Concentration | Lot dependent |
| Buffer/Preservatives | 0.01M PBS, pH 7.4. |
| Purity | >95% as determined by SDS-PAGE. |
| Purification | Protein A/G purified from cell culture supernatant. |
| Tested applications | IHC, WB |
| Antibody Type | Primary Antibody |
| Clone Number | BN349X4 |
| Endotoxins | Please contact with the lab for this information. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Expiration Date

12 months from date of receipt.



SDS-PAGE for Anti-GAD65/67 Antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)