

Product Datasheet

HCV NS1/gp68/gp70/Envelope glycoprotein E2 Human Monoclonal Antibody (orb2961020)

Catalog Number	orb2961020
Category	Antibodies
Description	Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody (BCR55) is a human monoclonal antibody detecting NS1 in ELISA, Neutralization.
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG1, kappa
Conjugation	Unconjugated
Reactivity	Virus
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
UniProt ID	Q5EG65, P26663, P27958
Tested applications	ELISA, NeA
Antibody Type	Recombinant Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	BCR55
Endotoxins	Please contact with the lab for this information.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2
Antibody binds with Envelope glycoprotein E2



Detects recombinant HCV glycoprotein E2 (Catalog No: orb2964080) in indirect ELISA.



SDS-PAGE for Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 Antibody.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Western blot analysis was performed using anti-HCV E2 monoclonal antibody at 1 µg/mL on various samples. Lane 1: recombinant HCV NS1/gp68/gp70/Envelope glycoprotein E2 (Catalog No: orb2964080), Lane 2: negative control.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)