

## Product Datasheet

### Recombinant Human FGL2 Protein, N-His (orb2962648)

|                             |                                                                                                                                                                           |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2962648                                                                                                                                                                |
| <b>Category</b>             | Proteins                                                                                                                                                                  |
| <b>Description</b>          | Recombinant Human FGL2 Protein, N-His                                                                                                                                     |
| <b>Form/Appearance</b>      | Lyophilized                                                                                                                                                               |
| <b>Buffer/Preservatives</b> | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| <b>Reconstitution</b>       | Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.                               |
| <b>Purity</b>               | >90% as determined by SDS-PAGE.                                                                                                                                           |
| <b>Protein Length</b>       | Asn24-Pro439                                                                                                                                                              |
| <b>MW</b>                   | 49.94 kDa                                                                                                                                                                 |
| <b>Expression System</b>    | E. coli                                                                                                                                                                   |
| <b>Endotoxins</b>           | Please contact with the lab for this information.                                                                                                                         |
| <b>Biological Origin</b>    | Homo sapiens (Human)                                                                                                                                                      |
| <b>Storage</b>              | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |
| <b>Note</b>                 | For research use only                                                                                                                                                     |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)