

Product Datasheet

Native Mouse IgG Protein (orb2966577)

Catalog Number	orb2966577
Category	Proteins
Description	Native Mouse IgG Protein
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Purity	>90% as determined by SDS-PAGE.
Protein Length	Full length.
MW	150 kDa
Expression System	Native
Endotoxins	Please contact with the lab for this information.
Biological Origin	Mus musculus (Mouse)
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558