

Product Datasheet

IgG Format Isotype Control for Bispecific Monoclonal Antibody (orb2970135)

Catalog Number	orb2970135
Category	Antibodies
Description	IgG Format Isotype Control for Bispecific Antibody is a human monoclonal antibody detecting Anti-Human HEL (HyHEL-10) & Anti-KLH Isotype Control in FCM, Isotype Control, In vivo Control . Highlights Multi-Application — Validated across multiple applications. Flow Cytometry — Validated for flow cytometry applications.
Clonality	Recombinant
Species/Host	Human
Isotype	IgG like (Tetravalent)
Conjugation	Unconjugated
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Tested applications	Control, FC
Antibody Type	Recombinant Antibody
Endotoxins	Please contact with the lab for this information.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)