

Product Datasheet

Mouse IgG2a, kappa Isotype Control Antibody (MOPC-173), PerCP (orb2970335)

Catalog Number	orb2970335
Category	Antibodies
Description	Mouse IgG2a, kappa Isotype Control Antibody (MOPC-173), PerCP
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a, kappa
Conjugation	PerCP
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
Purification	Affinity chromatography.
Tested applications	Control, FC
Endotoxins	Please contact with the lab for this information.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)