

Product Datasheet

Human CD85j/LILRB1/ILT2 Antibody (SAA0402), PE (orb2970568)

Catalog Number	orb2970568
Category	Antibodies
Description	Human CD85j/LILRB1/ILT2 Antibody (SAA0402), PE
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a, kappa
Conjugation	PE
Reactivity	Human
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
Tested applications	ELISA, FC
Endotoxins	Please contact with the lab for this information.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558