

Product Datasheet

CD184/CXCR4 Mouse Monoclonal Antibody (PE) (orb2970714)

Catalog Number	orb2970714
Category	Antibodies
Description	Anti-Human CD184/CXCR4 Antibody (BUCC2280), PE is a mouse monoclonal antibody detecting Stromal cell-derived factor 1 receptor in FCM. Suitable for Human. Highlights Flow Cytometry: Validated for flow cytometry applications.
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	PE
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
UniProt ID	P61073
Tested applications	FC
Clone Number	BUCC2280
Endotoxins	Please contact with the lab for this information.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Expiration Date

12 months from date of receipt.



Flow-cytometry using PE anti-human CD184 antibody. Jurkat cells were stained with an irrelevant antibody (Blue Histogram) or an PE anti-human CD184 monoclonal antibody (Catalog: orb2970714, Yellow Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, and cells analysed on a NovoCyte Flow Cytometer.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)