

Product Datasheet

Rat Monoclonal Antibody (PerCP) (orb2972240)

Catalog Number	orb2972240
Category	Antibodies
Description	Rat IgG2b, kappa Isotype Control Antibody (RTK4530), PerCP [RTK4530] is a rat monoclonal antibody in FCM, Isotype Control. Highlights Flow Cytometry — Validated for flow cytometry applications. Fluorescent Label — Direct conjugate for streamlined workflows.
Clonality	Monoclonal
Species/Host	Rat
Isotype	IgG2b, kappa
Conjugation	Unconjugated
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
Tested applications	Control, FC
Clone Number	RTK4530
Endotoxins	Please contact with the lab for this information.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558