

Product Datasheet

Human CD1a Antibody (2113), APC (orb2973190)

Catalog Number	orb2973190
Category	Antibodies
Description	Human CD1a Antibody (2113), APC
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1, kappa
Conjugation	APC
Reactivity	Human
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
Tested applications	FC
Endotoxins	Please contact with the lab for this information.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)