

Product Datasheet

CD25/IL2RA Rat Monoclonal Antibody (APC) (orb2973275)

Catalog Number	orb2973275
Category	Antibodies
Description	Anti-Mouse CD25/IL2RA Antibody (PC61/PC61.5.3), APC [PC61/PC61.5.3] is a rat monoclonal antibody detecting Interleukin-2 receptor subunit alpha in FCM . Suitable for Mouse . Highlights Flow Cytometry — Validated for flow cytometry applications. Fluorescent Label — Direct conjugate for streamlined workflows.
Clonality	Monoclonal
Species/Host	Rat
Isotype	IgG1, lambda
Conjugation	APC
Reactivity	Mouse
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
UniProt ID	P01590
Tested applications	FC
Clone Number	PC61/PC61.5.3
Endotoxins	Please contact with the lab for this information.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)