

Product Datasheet

Rat IgG2b Monoclonal Antibody (FITC) (orb2974255)

Catalog Number	orb2974255
Category	Antibodies
Description	Rat IgG2b Monoclonal Antibody (FITC)
Clonality	Monoclonal
Species/Host	Rat
Isotype	IgG2b, kappa
Conjugation	FITC
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	0.01M PBS, pH 7.4, 0.2% rAlbumin, 0.05% Proclin 300.
Tested applications	Control, FC
Clone Number	ES26-5E12.4
Endotoxins	Please contact with the lab for this information.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)