

Product Datasheet

CD41/ITGA2B/CD61/ITGB3 Mouse Monoclonal Antibody (orb2974557)

| | |
|-----------------------------|--|
| Catalog Number | orb2974557 |
| Category | Antibodies |
| Description | Anti-Human CD41/ITGA2B & CD61/ITGB3 Antibody (OPG2) [OPG2] is a mouse monoclonal antibody detecting CD41 in FCM. Suitable for Human. Highlights Flow Cytometry: Validated for flow cytometry applications. |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Isotype | IgG1, kappa |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Buffer/Preservatives | 0.01M PBS, pH 7.4, 0.09% Sodium azide. |
| Purity | >95% as determined by SDS-PAGE. |
| Purification | Protein A/G purified from cell culture supernatant. |
| UniProt ID | P05106, P08514 |
| Tested applications | FC |
| Clone Number | OPG2 |
| Endotoxins | Please contact with the lab for this information. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)