

Product Datasheet

CAT/Catalase (orb2975463)

Catalog Number	orb2975463
Category	Proteins
Description	CAT/Catalase expressed in E. coli . Purity: >90% by SDS-PAGE. Highlights E. coli Expression — High-yield, cost-effective production. High Purity — >90% purity verified by SDS-PAGE.
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS, pH 7.4.
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Purity	>90% as determined by SDS-PAGE.
Protein Length	CAT / Catalase is cloned from human cell lines and expressed in E.coli.
MW	62.07 kDa
Expression System	E. coli
Endotoxins	Please contact with the lab for this information.
Biological Origin	Human
Expression Region	CAT / Catalase is cloned from human cell lines and expressed in E.coli.
Storage	-20°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558