

Product Datasheet

Human ST2 ELISA Kit (orb2975521)

Catalog Number	orb2975521
Category	Assays and Kits
Description	Human ST2 ELISA Kit— Quantitative sandwich ELISA for Human ST2 ELISA Kit . Sample: Plasma, Serum . Range: 31.25 - 2, 000 pg/mL . Sensitivity: 21.22 pg/mL . Highlights Sandwich Format — High specificity, reliable quantification. High Sensitivity — Optimized for low-abundance target detection. Ready-to-Use — Complete kit with pre-coated plates included.
Tested applications	ELISA
Assay Type	Quantitative
Range	31.25 - 2,000 pg/mL
Sensitivity	21.22 pg/mL
Sample Types	Plasma, Serum
Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)