

## Product Datasheet

### Romiplostim ELISA Kit (orb2975554)

<b>Catalog Number</b>	orb2975554
<b>Category</b>	Assays and Kits
<b>Description</b>	<p>PRINCIPLE OF THE ASSAY This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for Romiplostim has been pre-coated onto a microplate. Standards or samples are pipetted into the wells and any Romiplostim present is bound by the immobilized antibody. After washing away any unbound substances, a biotin-labeled antibody specific for Romiplostim is added to the wells. After washing away any unbound substances, Streptavidin-HRP is added to the wells. Following a wash to remove any unbound enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of Romiplostim bound in the initial step. The color development is stopped and the intensity of the color is measured.</p>
<b>Assay Type</b>	Quantitative
<b>Range</b>	11.72 - 750 pg/mL
<b>Sensitivity</b>	11.31 pg/mL
<b>Sample Types</b>	Plasma, Serum
<b>Storage</b>	When the kit was stored at the recommended temperature for 12 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.
<b>Note</b>	For research use only

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558