

Product Datasheet

Anti-Secukinumab Neutralizing Antibody ELISA kit (orb2975636)

Catalog Number	orb2975636
Category	Assays and Kits
Description	Anti-Secukinumab Neutralizing Antibody ELISA kit— Quantitative sandwich ELISA for Anti-Secukinumab Neutralizing Antibody . Sample: Plasma, Serum . Range: 15.625-1, 000 ng/ml . Highlights Sandwich Format — High specificity, reliable quantification. Ready-to-Use — Complete kit with pre-coated plates included.
UniProt ID	Q16552
Tested applications	ELISA
Assay Type	Quantitative
Range	15.625-1,000 ng/ml
Sample Types	Plasma, Serum
Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)