

Product Datasheet

Bimagrumab ELISA Kit (orb2975740)

Catalog Number	orb2975740
Category	Assays and Kits
Description	Bimagrumab ELISA Kit: Quantitative sandwich ELISA for Bimagrumab. Sample: Plasma, Serum. Range: 125 - 4, 000 ng/mL. Sensitivity: 97.28 ng/mL. Highlights Sandwich Format: High specificity, reliable quantification. High Sensitivity: Optimized for low-abundance target detection. Ready-to-Use: Complete kit with pre-coated plates included.
UniProt ID	Q13705, P27037
Tested applications	ELISA
Assay Type	Quantitative
Range	125 - 4,000 ng/mL
Sensitivity	97.28 ng/mL
Sample Types	Plasma, Serum
Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)