

## Product Datasheet

### Sacituzumab ELISA Kit (orb2975790)

|                            |                                                                                                                                                                                                                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb2975790                                                                                                                                                                                                                                                                                                |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                                                                                                                                           |
| <b>Description</b>         | The Sacituzumab ELISA Kit is designed for the quantitative measurement of Sacituzumab in plasma and serum samples. The assay has a detection range of 78.13 - 5,000 ng/mL and a sensitivity of 40.3 ng/mL.                                                                                                |
| <b>UniProt ID</b>          | <b>P09758</b>                                                                                                                                                                                                                                                                                             |
| <b>Tested applications</b> | ELISA                                                                                                                                                                                                                                                                                                     |
| <b>Assay Type</b>          | Quantitative                                                                                                                                                                                                                                                                                              |
| <b>Range</b>               | 78.13 - 5,000 ng/mL                                                                                                                                                                                                                                                                                       |
| <b>Sensitivity</b>         | 40.3 ng/mL                                                                                                                                                                                                                                                                                                |
| <b>Sample Types</b>        | Plasma, Serum                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>             | When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                     |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                                                                                            |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)