

Product Datasheet

Labetuzumab ELISA Kit (orb2975802)

Catalog Number	orb2975802
Category	Assays and Kits
Description	The Labetuzumab ELISA Kit is designed for the quantitative measurement of Labetuzumab in plasma and serum samples. The assay has a detection range of 156.25 - 10,000 ng/mL and a sensitivity of 53.67 ng/mL.
UniProt ID	P06731
Tested applications	ELISA
Assay Type	Quantitative
Range	156.25 - 10,000 ng/mL
Sensitivity	53.67 ng/mL
Sample Types	Plasma, Serum
Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)