

Product Datasheet

Besufetamig Human Monoclonal Antibody (orb2975859)

Catalog Number	orb2975859
Category	Antibodies
Description	besufetamig is a research-grade recombinant antibody targeting CD274. Highlights Research Grade: For PK/PD studies, assay development, and ADA research.
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG1, Kappa
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Purity	>95% as determined by SDS-PAGE.
UniProt ID	P07766, Q9NZQ7
Tested applications	ELISA, FA, FACS, In vivo
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.



SDS-PAGE for Research Grade besufetamig.



Western blot analysis was performed using anti-CD3E monoclonal antibody at 1 $\mu\text{g}/\text{mL}$ on various samples. Lane 1: recombinant human CD3E.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558