

Product Datasheet

Pasritamig Human Monoclonal Antibody (orb2975932)

Catalog Number	orb2975932
Category	Antibodies
Description	Pasritamig Human Monoclonal Antibody
Clonality	Monoclonal
Species/Host	Human
Isotype	Ig(H-gamma1_L-kappa)_scFvkh-G1(h-CH2-CH3)
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Purity	>95% as determined by SDS-PAGE.
Tested applications	ELISA, FA, FACS, In vivo
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558