

Product Datasheet

Velinotamig Human Monoclonal Antibody (orb2975961)

Catalog Number	orb2975961
Category	Antibodies
Description	Velinotamig Human Monoclonal Antibody is an unconjugated monoclonal antibody. It shows reactivity with Human samples. It is suitable for use in ELISA, FA, FACS, In vivo applications. The antibody is supplied in liquid form. It is Purity is >95% as determined by SDS-PAGE.
Target	CD3E & CD269/TNFRSF17/BCMA
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG1, Kappa
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	Lot dependent
Reconstitution	Reconstitute in sterile water for a stock solution. A copy of datasheet will be provided with the products, please refer to it for details.
Purity	>95% as determined by SDS-PAGE.
UniProt ID	P07766, Q02223
Tested applications	ELISA, FA, FACS, In vivo

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Biosimilar Antibody
CAS Number	2848586-49-0
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



SDS-PAGE for Research Grade velinotamig.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558