

Product Datasheet

Recombinant human CSF3R protein, C-His (HEK293) (orb2978568)

Catalog Number	orb2978568
Category	Proteins
Description	Recombinant human GCSF Receptor protein, C-His (HEK293), with a concentration of >1mg/ml, expressed in HEK293 cell, with a molecular weight of 68.4 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	Granulocyte colony-stimulating factor receptor
Tag	C-His
Form/Appearance	Liquid
Concentration	>1mg/ml
Buffer/Preservatives	PBS (pH=7.4)
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	>90% as determined by SDS-PAGE
Purification	AC
Protein Sequence	25-621/836
UniProt ID	Q99062
MW	68.4 kDa
Expression System	HEK293 cell

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)