

Product Datasheet

Anti-CD3 & BCMA Reference Antibody (teclistamab Biosimilar) (orb2980125)

Catalog Number	orb2980125
Category	Antibodies
Description	Anti-CD3 & BCMA Reference Antibody (teclistamab Biosimilar)
Clonality	Monoclonal
Isotype	IgG-like
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Lyophilized
Buffer/Preservatives	100 mM Proline/Acetic acid 20mM Arginine pH 5.0
Reconstitution	Dissolve with sterile ddH2O
Purification	Protein A
MW	143.64 kDa
Tested applications	ELISA, FA, FACS, In vivo
Antibody Type	Primary Antibody
Expression System	CHO cell
Endotoxins	0.001 EU/ug

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)