

Product Datasheet

Mouse Protein Kinase, AMP Activated Alpha 1 (PRKAa1) ELISA Kit (orb2980242)

Catalog Number	orb2980242
Category	Assays and Kits
Description	The Mouse Protein Kinase, AMP Activated Alpha 1 (PRKAa1) ELISA Kit features a dynamic range of 78.1-5000pg/ml and is suitable for Serum, plasma, tissue homogenates or other biological fluids.
Reactivity	Mouse
Range	78.1-5000pg/ml
Sample Types	Serum, plasma, tissue homogenates or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)