

Product Datasheet

Human 14-3-3 Protein Epsilon (YWHAE) ELISA Kit (orb2980251)

Catalog Number	orb2980251
Category	Assays and Kits
Description	The Human 14-3-3 Protein Epsilon (YWHAE) ELISA Kit features a dynamic range of 78.1-5000pg/mL and is suitable for Serum, plasma, tissue homogenates or other biological fluids.
Reactivity	Human
Range	78.1-5000pg/mL
Sample Types	Serum, plasma, tissue homogenates or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)