

Product Datasheet

Rat S100A1 Calcium Binding Protein (S100A1) ELISA Kit (orb2981106)

Catalog Number	orb2981106
Category	Assays and Kits
Description	The Rat S100A1 Calcium Binding Protein (S100A1) ELISA Kit features a dynamic range of 1.56-100ng/mL and is suitable for Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Reactivity	Rat
Range	1.56-100ng/mL
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)