

Product Datasheet

Rat Platelet Derived Growth Factor Receptor Alpha (PDGFRa) ELISA Kit (orb2981567)

Catalog Number	orb2981567
Category	Assays and Kits
Description	The Rat Platelet Derived Growth Factor Receptor Alpha (PDGFRa) ELISA Kit features a dynamic range of 0.156-10ng/mL and is suitable for Tissue homogenates, cell lysates or other biological fluids.
Reactivity	Rat
Range	0.156-10ng/mL
Sample Types	Tissue homogenates, cell lysates or other biological fluids.
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)