

## Product Datasheet

# Mouse Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit (orb2983089)

<b>Catalog Number</b>	orb2983089
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit features a dynamic range of 0.156-10ng/mL and is suitable for Tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Sample Types</b>	Tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)