

## Product Datasheet

# ULK1 (Phospho-S556) Rabbit Monoclonal Antibody (orb2988906)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb2988906  |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | The Recombinant ULK1 (Phospho-S556) Rabbit mAb Antibody is suitable for WB. It is a Monoclonal, Unconjugated antibody which raised against KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S556 of human ULK1 protein. The exact sequence is proprietary. Purification: The antibody was purified by immunogen affinity chromatography. |
| <b>Target</b>               | ULK1  |
| <b>Clonality</b>            | Monoclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Reactivity</b>           | Human, Mouse  |
| <b>Form/Appearance</b>      | Liquid  |
| <b>Buffer/Preservatives</b> | PBS, pH 7.4, containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.   |
| <b>Purification</b>         | The antibody was purified by immunogen affinity chromatography.   |
| <b>Immunogen</b>            | KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S556 of human ULK1 protein. The exact sequence is proprietary.  |
| <b>UniProt ID</b>           | <b>070405, 075385</b>   |
| <b>Tested applications</b>  | WB  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |   |
|------------------------|---|
| <b>Dilution range</b>  | WB (1/500 - 1/1000)   |
| <b>Specificity</b>     | Recognizes endogenous levels of ULK1 protein only when phosphorylated at S556   |
| <b>Antibody Type</b>   | Primary Antibody  |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Entrez</b>          | <b>8408, 22241</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.   |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)