

Product Datasheet

Recombinant FGFR2 Rabbit mAb Antibody (orb2989747)

Catalog Number	orb2989747
Category	Antibodies
Description	The Recombinant FGFR2 Rabbit mAb Antibody is suitable for IHC. It is a Monoclonal, Unconjugated antibody which raised against Recombinant fusion protein of human FGFR2. The exact sequence is proprietary. Purification: The antibody was purified by immunogen affinity chromatography.
Clonality	Monoclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid in PBS, pH 7.3, 50% glycerol, and 0.05% Proclin300.
Purification	The antibody was purified by immunogen affinity chromatography.
Immunogen	Recombinant fusion protein of human FGFR2. The exact sequence is proprietary.
UniProt ID	P21802
Tested applications	IHC
Dilution range	IH (1/100 - 1/300)
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Explore. Bioreagents.

Biorbyt.com

Entrez

2263

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44%280%201223%20859-353) | Fax: [+1 \(415\) 651-8558](tel:+1%28415%20651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1%28415%20906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1%28415%20651-8558)