

Product Datasheet

Human Siglec-2/CD22 ELISA Kit (orb2991724)

Catalog Number	orb2991724
Category	Assays and Kits
Description	Human Siglec-2/CD22 ELISA Kit (96 Tests). Quantitate Human Siglec-2/CD22 in cell culture supernatants, serum, plasma(heparin, EDTA). Sensitivity: 15pg/ml. The brand indicates this is a premium quality ELISA kit. Each kit delivers precise quantification, high sensitivity, and excellent reproducibility. Only our most reliable and effective kits qualify as , guaranteeing top-tier results for your assays.
Reactivity	Human
UniProt ID	P20273
MW	95348 Da
Tested applications	ELISA
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	93.8-6,000 pg/mL
Sensitivity	15 pg/mL
Sample Types	cell culture supernatants, serum, plasma(heparin, EDTA)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)